
CRISTAVISC c®

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
DESCRIPTION

CRISTAVISC c® Ophthalmic Viscosurgical Device (OVD) is a viscoelastic sodium hyaluronate (155 % gel of non-animal origin, with a high molecular weight (2 to 3 MDa) and incorporating an antioxidant (mannitol). It is transparent, odorless, sterile and non-pyrogenic, with physiological pH and osmolality. The gel is packaged in a pre-filled syringe and sterilised using autoclave. The blister is sterilised with ethylene oxide. This OVD is for single use. Each box contains one 1ml **CRISTAVISC c®** pre-filled syringe, one single-use sterile cannula, an instruction leaflet and labels with the batch number.

CRISTAVISC c® has an osmolality of 250-400 mOsmol/L, a pH of 6.8 to 7.6 and a viscosity of approximately 170 Pas at a shear rate of 0.01 s⁻¹ (T°=25°C).

INDICATIONS

CRISTAVISC c®	An Ophthalmic Viscosurgical Device indicated for use as a surgical aid in ophthalmic anterior segment surgery (protects the corneal endothelium and maintains the intraocular space).
---------------	---

MODE OF ACTION

CRISTAVISC c® is easy to install and to remove from the anterior chamber at the end of the operation. During the injection, the viscosity decreases due to the shearing in the cannula. The original high viscosity is then restored after installation. Due to its specific viscoelastic properties, **CRISTAVISC c®** maintains a deep anterior chamber and protects the corneal endothelium throughout the surgery. The viscoelastic and transparency properties of the product enable the atraumatic and precise handling of eye tissues during eye surgery. The mannitol in **CRISTAVISC c®** acts as a free radical scavenger, particularly for hydroxyls, thus decreasing oxidative stress during eye surgery. Its presence limits the destruction of the sodium hyaluronate and, as a result, helps it to retain its rheological properties especially during phacemulsification. Once the surgery is finished, **CRISTAVISC c®** must be completely removed by irrigation and/or aspiration. Any potential residues will be flushed out of the anterior chamber naturally via the trabecular meshwork and Schlemm's canal.

CONTRAINDICATIONS

The OVD must not be used:
- In patients with a known hypersensitivity to one of the ingredients.

INTERACTIONS

CRISTAVISC c® must not come into contact with ophthalmic preparations containing quaternary ammonium compounds as a preservative or with surgical instruments rinsed with this type of solution. There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds, such as benzalkonium chloride solutions, which can cause the formation of a precipitate.

PRECAUTIONS FOR USE

- The product must be used for intraocular use only.
- Patients with pre-existing glaucoma, severe myopia, diabetic retinopathy or uveitis prior to surgery are at higher risk of increased intraocular pressure.
- The OVD has not been evaluated in pregnant or breastfeeding women or children.
- Only the cannulas provided with the OVD should be used for the injection, as the combination of these two devices has been validated.
- Do not use if the packaging has been damaged (syringe, blister pack, lid).
- Inject immediately after opening.
- At the end of the surgical procedure, it is necessary to completely remove the OVD from the eye by thorough irrigation and/or aspiration. However there is an occasional risk of blockage of the trabecular meshwork, causing a temporary post-operative rise in intraocular pressure.
- Discard the syringe and the residual product when the injection is complete, and dispose of the cannula in an appropriate container, in accordance with the current national regulations.
- The product is for single use. Do not reuse.
- A product must be used only for one single patient to prevent any risk of cross-contamination.
- After opening, the product can never be re-sterilised, even if the OVD has not been injected.

METHOD OF ADMINISTRATION

- The treatment should be carried out in an operating room by a legally qualified practitioner with good knowledge of the ophthalmic anterior segment surgery.
- Before carrying out the surgery, it is essential:
- to thoroughly confirm the patient's medical history with her/him;
- to explain the contraindications and potential adverse effects related to the procedure with the patient.

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CRISTAVISC c®	Natriumhyaluronat Fosfatbuffer og mannitol - pH 72 1 ml forfylt sprøjte	155 mg q.s 1 g
---------------	---	-------------------

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
- Before injection:
- Verify the integrity of the cannula;
- Verify that the gel is not cloudy;
- Observe the aspect and handling rules related to this type of intervention.

Assembly-instructions

1. Remove the plastic end cap of the syringe by twisting it off carefully.
 2. Hold the ring at the end of the syringe and firmly insert the hub of the cannula into the Luer lock end of the syringe, by screwing until you feel some strong counter pressure.
- Avoid trapping air bubble.
- Slowly inject, using the cannula provided.
- During the procedure, additional OVD can be installed to maintain the anterior chamber or replace any viscoelastic lost during the surgery.
- When injected in the anterior chamber of the eye, the amount of the OVD should be adjusted according to the volume of the aqueous humour. Do not overfill the eye chamber. Overfilling the anterior segment of the eye with the OVD may cause increased intraocular pressure, glaucoma, or other ocular damage.
- If the cannula is blocked, do not increase the pressure required for injection. Replace the cannula.
- At the end of the surgical procedure, completely remove the OVD from the eye by thorough irrigation and/or aspiration.
- Carefully monitor the intraocular pressure, especially during the immediate post-operative period. If significant post-operative intraocular pressure rises are observed, appropriate corrective therapy must be initiated.

ADVERSE EFFECTS

The practitioner must inform the patient of any potential adverse effects that may occur following the use of the OVD, in particular (this list is not exhaustive):
- Post-operative transient increase in intraocular pressure if the OVD is not removed as completely as possible. A careful monitoring must be set up to manage any increased post-operative intraocular pressure and to reduce the likelihood of occurrence of secondary glaucoma or other ocular damage.
- Rarely post-operative inflammatory reactions such as toxic anterior segment syndrome (TASS), endophthalmitis, as well as corneal oedema have been reported in literature following injections of sodium hyaluronate. Nevertheless, no relationship with the product has been established.
- Rare cases of hypersensitivity have been reported in literature following injections of sodium hyaluronate. The patient must be informed about this.
- The patient must inform the practitioner of any adverse effects mentioned above, or the appearance of any other adverse effect. The practitioner must report it to the reseller or manufacturer as soon as possible and should carry out an appropriate care.

STORAGE - EXPIRY DATE

- The OVD must be used before the expiration date indicated on the packaging and on the syringe.
- The OVD must be stored in its original packaging, at a temperature between 2°C and 25°C [36°F and 77°F], away from frost and light.
- If storage conditions are not observed, the OVD's performance can be compromised.

Caution. Consult the instruction leaflet.	Temperature limits. Store between 2°C and 25°C [36°F and 77°F].
Expiry date. Use before the date shown.	Keep away from sunlight.
Batch code.	Fragile, handle with care.
LOT	Do not reuse.
INTÉGRITÉ	Do not use if the packaging has been damaged (syringe, blister pack, lid).
STÉRILISÉ	Manufacturer.
CE	The CE marking was obtained in compliance with Directive 93/42/EEC relative to medical devices.
0344	CRISTAVISC c® obtained CE marking in 2010. 0344 corresponds to the number of the Notified Body.

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