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

DESCRIZIONE
PHYSIOVISC® INTEGRAL Ophthalmic Viscosurgical Device (OVD) is a viscoelastic sodium hyaluronate (1,55 % of non-animal origin, with a high molecular weight (2 to 3 MDa) and incorporating an antioxidant (mannitol). It is transparent, odourless, sterile and non-pyrogenic, with physiological pH and compatibility.
The gel is packaged in a pre-filled syringe and sterilised using autoclave. The blister is sterilised with ethylene oxide. This OVD is single-use.
Each box contains one **PHYSIOVISC® INTEGRAL** pre-filled syringe, one single-use sterile cannula, an instruction leaflet and labels with the batch number.

PHYSIOVISC® INTEGRAL has a compatibility of 250-400 mOsmol/L a pH of 6.8 to 7.6 and a viscosity of approximately 170 Pa.s at a shear rate of 0.01s⁻¹ (T⁺=25 °C).

INDICAZIONI

An Ophthalmic Viscosurgical Device indicated for use as a surgical aid in ophthalmic anterior segment surgery (protects the corneal endothelium and maintains the intraocular space).

MODE D'ACTION

PHYSIOVISC® INTEGRAL is easy to install and to remove from the anterior chamber at the end of the operation. During the injection, the viscosity decreases due to the shearing in the cannula. The original high viscosity is then restored after installation.
Due to its specific viscoelastic properties, **PHYSIOVISC® INTEGRAL** maintains a deep anterior chamber and protects the corneal endothelium throughout the surgery. The viscoelastic and transparency properties of the product ensure the atraumatic and precise handling of eye tissues during eye surgery.
The mannitol in **PHYSIOVISC® INTEGRAL** acts as a free radical scavenger, particularly for hydroxyl radicals, thus decreasing oxidative stress during eye surgery. Its presence limits the destruction of the sodium hyaluronate and, as a result, helps it to retain its rheological properties especially during pharmacemulsification.
Once the surgery is finished, **PHYSIOVISC® INTEGRAL** must be completely removed by irrigation and/or aspiration. Any potential residues will be flushed out of the anterior chamber naturally via the trabecular meshwork and Schlemm's canal.

CONTRAINDICATIONS

- Do not use if the patient has a known hypersensitivity to one of the ingredients.

INTERACTIONS

PHYSIOVISC® INTEGRAL must not come into contact with ophthalmic preparations containing quaternary ammonium compounds as a preservative or with surgical instruments rinsed with this type of solution. There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds, such as benzalkonium chloride solutions, which can cause the formation of a precipitate.

PRECAUTIONS FOR USE

- The product must be used for intraocular use only.
- Patients with pre-existing glaucoma, severe myopia, diabetic retinopathy or uveitis prior to surgery are at a higher risk of increased intraocular pressure.
- The OVD has not been evaluated in pregnant or breastfeeding women or children.
- Only the cannulas provided with the OVD should be used for the injection, as the combination of these two devices has been validated.
- Do not use if the packaging has been damaged (syringe, blister pack, lid).
- Inject immediately after opening.
- At the end of the surgical procedure, it is necessary to completely remove the OVD from the eye by thorough irrigation and/or aspiration; otherwise there is an occasional risk of blockage of the trabecular meshwork, causing a temporary post-operative rise in intraocular pressure.
- Discard the syringe and the residual product when the injection is complete, and dispose of the cannulas in an appropriate container, in accordance with the current national regulations.
- The product is for single use. Do not reuse.
- A product must be used only for one single patient to prevent any risk of cross-contamination.
- After opening, the product can never be re-sterilised, even if the OVD has not been injected.

METHOD OF ADMINISTRATION

- The treatment should be carried out in an operating room by a legally qualified practitioner with good knowledge of the ophthalmic anterior segment surgery.
- Before carrying out the surgery, it is essential:
- to thoroughly confirm the patient's medical history with her/him;
- to explain the contraindications and potential adverse effects related to the procedure with the patient.

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